









REQUIREMENTS FOR MICROVISION SOFTWARE INSTALLATION

The Microvision Classic and Advance software is compatible with all c-arms that have a 525/625 line and 50Hz/60Hz composite video signal.

To avoid short-term hardware failures, we recommend using completely new and unused hardware.

Monitor	
	Led TV - 32" Led Monitor compatible with HD resolution 1366x768p
Computer	
	Intel processors From 2nd to 7th generation, Intel Core I5 From 8th generation onwards, Intel Core I3 From 10th generation onwards, Intel Pentium Gold G6400
	AMD processors AMD Ryzen 3 3200 en onwards
	Ram memory 4 GB onwards
	diskdrive SSD 240GB onwards
	O.S. Windows 10 64 bit

	<p>AVerMedia PCIe capture card C725 or CE310B Sample Link 1 Sample Link 2</p>
	<p>EZCAP170 or EZCAP172 USB video grabber Sample Link 1</p>
	<p>"Easycap" or similar USB video grabber with UTV007 chipset</p>

INSTALLING MICROVISION SOFTWARE

Our goal is for you to achieve the best experience and obtain the best results in the testing process, which is why the software installation will be carried out exclusively by us remotely through Teamviewer or Anydesk.

From our place we will provide you with full support to make the necessary adjustments and configurations in the adaptation of the Microvision Software with its C-Arm, for this reason I kindly request that you contact us via Whatsapp or Email contacto@microvision-software.com and we will work together the day and time agreed upon.

CONNECTION EXAMPLE FOR TESTING WITH SIREMOBIL COMPACT

Step 1: Open the monitor trolley and detect the video signal from the camera



Step 2: Use a cable to connect the video from the camera to the USB grabber



Step 3: Connect the video cable to the USB grabber and plug it into the PC



Step 4: Contact me via TeamViewer and Whatsapp to configure the software.

